






Mailing Kit Instructions for Forensic Toxicology

Components:

1. Specimen Collection Instructions and Testing Requisition.
2. Two (2) purple-top, (EDTA) or two (2) grey-top (NaF) blood specimen tubes.
3. Absorbent tube sleeve.
4. Zipper lock biohazard bag.
5. Corrugated shipping carton.
6. Hazardous Material Shipping Label (UN3373).

Not Supplied:

1. Dry ice.
2. Sealing/Packing tape.

Step 1		<ol style="list-style-type: none"> a. Place a tube in one slot of the absorbent sleeve. b. Repeat for the 2nd and 3rd blood tubes (if needed).
Step 2		<ol style="list-style-type: none"> a. Place absorbent sleeve with tubes in the biohazard bag. b. Seal the bag.
Step 3		<ol style="list-style-type: none"> a. Place the specimen bag and testing requisition in the corrugated shipping carton.
Step 4		<ol style="list-style-type: none"> a. Close the carton end flaps and insert the carton tab. b. Place tape over the flap to seal the box.
Step 5		<ol style="list-style-type: none"> a. Place the UN3373 label on outside of a biological shipping container. (if shipping via US postal service). b. Seal the outer box along all seams for delivery via preferred carrier (i.e. UPS, USPS, or FedEx).

***Certification that this Product satisfies the requirements of applicable Domestic and International regulations for the shipment of diagnostic specimens is contingent upon the proper use and assembly of all supplied components. Under applicable regulations, the shipper is responsible for selecting proper packaging and final closure. Shipper should read the applicable regulations and make its own determination as to what are the requirements for packaging a particular shipment.

For questions, comments, or concerns contact MolecularDx at 814-361-6997/info@molecdx.com